

# Oakworks CFPM100 Imaging Table w/Rectangular Top

## FACTORY NEW

The Oakworks CFPM100 Imaging Table is a single-tower lift system designed for pain management practices, surgical centers, and multi-disciplinary medical facilities. It features a 500-pound weight capacity and is compatible with all mobile C-Arm imaging systems, providing flexibility for various facility configurations.



Stock Photo - Options May Vary

The table offers a height range of 22 inches to 40 inches, accommodating ergonomic working requirements while enabling safe wheelchair transfers at its lowest setting. The rectangular carbon fiber top measures 22 by 84 inches and provides 56 inches of metal-free imaging space for fluoroscopy-guided procedures.

The CFPM100 incorporates a soft-touch table base material that protects against C-Arm damage and scrapes. It features a state-of-the-art foot control that can remain in the base or be removed and repositioned to either side of the table for operator convenience.

## FEATURES

- Rectangular Top
- A single-tower lift system
- Electric Height Adjustment with a range of 22" to 40"
- Weight capacity: 500 lbs distributed evenly
- Compatible with all mobile C-Arm imaging systems
- Soft-touch table base material
- Durable base material Protects and Minimizes c-arm scrapes and damage
- Included with Table: 22" x 84" tabletop pad along, with and Safety Strap

## SPECIFICATIONS

- Carbon Fiber Rectangular Top Dimensions: 22"x 84"
- Imaging Space: 56" of Metal-Free Space
- Patient Weight Capacity: 500 lbs distributed evenly
- Electronics: 110 Volt Motor
- 4" Locking Casters: Locking / Unlocking casters for mobile or permanent positioning
- Table Weight: 343 lbs
- Frame: Single Tower Steel Base
- Warranty: 3 years parts, 2 years labor

*\*Brochure and images are a general overview of the product. Actual product options may vary.*