

Oakworks CFPM200 Imaging Table with Rectangular Top

FACTORY NEW

The Oakworks CFPM200 Imaging Table with Rectangular Top is a high-performance 2-motion table designed for pain management and surgical imaging environments. The carbon fiber rectangular top delivers exceptional radiolucency for clear and precise imaging results. Its removable Comfort Foam® pad ensures optimal patient comfort without compromising image quality, making it ideal for extended procedures.



Stock Photo - Options May Vary

Engineered for stability and ease of use, the CFPM200 features a wide height range that enhances safety during patient transfers and allows clinicians to work at ergonomically favorable positions. The electronic lift tower enables smooth vertical motion, while the robust lifting column provides rock-solid support during demanding procedures. An emergency stop switch offers an added layer of operational safety by cutting power immediately when necessary.

The table is designed for both flexibility and secure positioning in clinical settings. Locking casters allow for easy mobility or stable, fixed placement as needed. All functions are controlled via hand-held remote, enabling quick and efficient adjustments. Two patient restraint straps secure individuals comfortably throughout the imaging process.

FEATURES

- Carbon fiber rectangular top with removable radiolucent Comfort Foam® table top pad
- All functions hand control
- Electronic lift tower
- Emergency stop switch to cut power from the table in an emergency
- Locking casters for mobile or permanent positioning
- (2) Patient safety straps
- Motion 1 - Electric Height Range 25"- 43" (64-109cm) Height Range without pad
- Motion 2 - Electric Lateral tilt $\pm 15^\circ$

SPECIFICATIONS

- Carbon Fiber Top - 22"x84" (56cmx214cm)
- 52" (133cm) Metal free imaging space
- Motion 1 - Electric Height Range 25"- 43" (64-109cm) Height Range without pad
- Motion 2 - Electric Lateral tilt $\pm 15^\circ$
- Weight Capacity - 500 lb (227kg)
- Product Weight - 475 lb (216kg)
- Locking Casters - 4in (10cm)
- 2 in. (5cm) Radiolucent Comfort Foam® table top pad

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com