

Oakworks CFPM201 Imaging Table with Integrated Headrest



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPM201 is a 2-motion C-Arm imaging table designed for pain management procedures. It features a carbon fiber integrated headrest top with removable radiolucent Comfort Foam® pad and a crescent face rest pad for prone positioning. The large, rigid lifting column provides exceptional stability during imaging.

The table offers a wide electric height range from 25" to 42" to support safe patient transfers and comfortable clinician positioning. Electric Trendelenburg tilt capability of +15°/-11° enables precise patient positioning for complex procedures. All motion functions are operated via hand control for efficient workflow.

Safety features include an emergency stop switch, locking casters for secure or mobile positioning, and two patient safety straps. The carbon fiber top measures 22"x84" with 56" of metal-free imaging space, supporting a 500 lb weight capacity.

FEATURES

- Carbon fiber integrated headrest top with removable radiolucent Comfort Foam® table top pad
- All functions hand control
- Electronic lift tower
- Emergency stop switch to cut power from the table in an emergency
- Locking casters for mobile or permanent positioning
- (2) Patient safety straps
- Crescent face rest pad (included with integrated headrest top)
- Motion 2 - Electric Trendelenburg tilt +15°/-11°

SPECIFICATIONS

- Carbon Fiber Top - 22"x84" (56cmx214cm)
- 56" (133cm) Metal free imaging space
- Motion 1 - Electric 25"- 42" (64-107cm) Height Range without pad
- Motion 2 - Electric Trendelenburg tilt +15°/-11°
- Weight Capacity - 500 lb (227kg)
- Product Weight - 475 lb (216kg)
- 4 in. (10cm) locking casters
- Warranty - 4 years parts, 3 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*