

# Oakworks CFPM201 Imaging Table with Rectangular Top



Stock Photo - Options May Vary

## FACTORY NEW

The Oakworks CFPM201 Imaging Table is a 2-motion table designed for medical imaging and pain management procedures. It features a carbon fiber top providing 56 inches of metal-free imaging space ideal for C-Arm use, combined with a 2-inch radiolucent Comfort Foam® pad for patient comfort. The table includes 4-inch locking casters and two safety straps for secure positioning and mobility.

The table offers a wide height range of 25 to 42 inches, enabling safe patient transfers and ergonomic working positions for clinicians. A rigid lifting column delivers outstanding stability during imaging procedures. Hand-held remote control manages Trendelenburg tilt adjustment from +15° to -11° without disrupting workflow.

With a 500 lb weight capacity and premium PVC-free TerraTouch™ fabric, the CFPM201 supports durable use in demanding clinical settings. The table meets CE marking, cETLus approval, CB certification, and FDA listing standards, backed by a 4-year parts and 3-year labor warranty.

## FEATURES

- All functions hand control
- Electronic lift tower
- Emergency stop switch to cut power from the table in an emergency
- Locking casters for mobile or permanent positioning
- (2) Patient safety straps
- Large height range table 25" - 42"
- 56" (133 cm) of metal free imaging space
- 500 lb patient weight capacity

## SPECIFICATIONS

- Carbon Fiber Top - 22"x84" (56cmx214cm)
- 56" (133cm) Metal free imaging space
- Motion 1 - Electric Height Range 25"- 42" (64-107cm) Height Range without pad
- Motion 2 - Electric Trendelenburg tilt +15°/-11°
- Weight Capacity - 500 lb (227kg)
- Product Weight - 475 lb (216kg)
- 2 in. (5cm) Radiolucent Comfort Foam® table top pad
- 4 in. (10cm) locking casters

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)