

Oakworks CFPM300 Imaging Table w/Integrated Headrest



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPM300 Imaging Table features three powered motions for height adjustment, Trendelenburg tilt, and lateral tilt, accommodating patients up to 500 lbs. The table offers an industry-leading adjustable height range of 18 inches (26" to 43") and includes 56 inches of metal-free imaging area for optimal visualization during procedures.

The table comes equipped with an integrated headrest, hand control, patient safety strap with Velcro closure, 4-inch casters, and a 1-inch comfort foam pad. The rectangular carbon fiber top measures 22 inches wide by 84 inches long, providing a stable platform for imaging and treatment applications.

Manufactured in the USA, the CFPM300 is FDA listed, cETL approved, and CE marked. It includes a 3-year parts warranty and 2-year labor warranty.

FEATURES

- Integrated Headrest
- 3 Powered Motions: Height, Trendelenburg Tilt, and Lateral Tilt
- Industry's Widest Adjustable Height Range of 18" (26" to 44")
- Patient Weight Capacity: 500 lbs (227 kg).
- 56" of Metal-Free Imaging Area
- Patient Safety Strap with Velcro® Closure
- Hand Control Included
- Made in the USA

SPECIFICATIONS

- Patient Weight Capacity: 500 lbs
- Rectangular Carbon Fiber Top: 22" (56 cm) wide x 84" (214 cm) long.
- Motion 1: Electric Height Range 26 in - 43 in (66-109 cm) Height Range without pad.
- Motion 2: Electric Trendelenburg Tilt +15°/-11°.
- Motion 3: Electric Lateral Tilt ±15°.
- Control: 1 Hand (3 Motions)
- Table Weight: 475 lbs (216 kg).
- Warranty: 3 years parts, 2 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)