

Oakworks CFPM301 Imaging Table: Rectangular Top



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPM301 is a rectangular imaging table designed for diagnostic and interventional procedures. It features a 24" x 80" carbon fiber top with 56 inches of metal-free imaging space, accommodating a 500 lb patient weight capacity. The table includes three powered motions: height adjustment from 26" to 44", Trendelenburg tilt of +15°/-11°, and 10" of longitudinal travel.

The powered height range enables seamless transitions from low-height wheelchair transfers to maximum elevation for C-Arm lateral imaging requirements. All motions are controlled via a single multi-function hand control for intuitive operation. The table comes standard with a 1-inch radiolucent comfort foam pad, patient safety straps, and 4-inch locking casters.

The CFPM301 meets stringent safety standards with FDA listing, CE marking, cETLus approval, and CB certification. It includes a comprehensive warranty covering parts and labor.

FEATURES

- Carbon Fiber Top - 56 in. Metal-free imaging space
- Top Dimensions - 24" Wide x 80" Long
- Motion 1 Electric Height Range 26" - 44" (without pad)
- Motion 2 Electric Trendelenburg tilt +15°/-11°
- Motion 3 Electric Longitudinal travel 10 in.
- Weight Capacity - 500 lb
- Included Features - 1 in. Radiolucent Comfort Foam™ tabletop pad
- Includes Patient safety straps (2)

SPECIFICATIONS

- Patient Weight Capacity: 500 lbs.
- Carbon Fiber Top: 24" W x 80" L
- Motion 1: height range 26" - 44" (without pad)
- Motion 2: Trendelenburg tilt +15/-11°
- Motion 3: 10" longitudinal travel
- Control: 1 hand (3 motion)
- Table Weight: 475 lbs.
- Warranty: 3 years parts, 2 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)