

# Oakworks CFPM302 Imaging Table with Integrated Headrest

## FACTORY NEW

The Oakworks CFPM302 Imaging Table with Integrated Headrest is designed for interventional procedures with a carbon fiber integrated headrest top that ensures optimal radiolucency and imaging access. The table includes a removable Comfort Foam® pad for patient support and comfort during lengthy procedures. Its 3-motor system allows smooth positioning for head, height, and Trendelenburg adjustments, with a wide height range that promotes safer patient transfers and ergonomic working conditions for clinicians.



*Stock Photo - Options May Vary*

The table features a sturdy electronic lift tower for enhanced stability and locking casters that allow easy mobility or fixed placement. An integrated face rest with crescent face pad adds convenience for prone positioning during imaging or treatments. Hand-controlled adjustments and an emergency stop switch provide enhanced safety and control, with two patient safety straps to secure the patient during procedures.

## FEATURES

- Carbon fiber integrated headrest top with removable radiolucent Comfort Foam® table top pad
- Integrated face rest with crescent face pad
- All functions hand control
- Electronic lift tower
- Emergency stop switch to cut power from the table in an emergency
- Locking casters for mobile or permanent positioning
- (2) Patient safety straps
- Motion 1 - Electric 25"- 43" (64-109cm) Height Range without pad

## SPECIFICATIONS

- Carbon Fiber Top - 22"x84" (56cmx214cm)
- 56" (133cm) Metal free imaging space
- Weight Capacity - 500 lb (227kg)
- Product Weight - 475 lb (216kg)
- 2 in. (5cm) Radiolucent Comfort Foam® table top pad
- 4 in. (10cm) locking casters
- Warranty - 4 years parts, 3 years labor
- Safety Listings - CE marked, cETLus approved, CB Certificate, FDA listed

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)