

# Oakworks CFPM301 Imaging Table with Integrated Headset



Stock Photo - Options May Vary

## FACTORY NEW

The Oakworks CFPM301 is a carbon fiber imaging table designed for advanced medical imaging procedures. It features a 22-inch wide by 84-inch long carbon fiber top with 56 inches of metal-free imaging space, enabling clear visualization during scans. The table supports up to 500 pounds and weighs 475 pounds.

Electric-powered adjustments provide versatile positioning capabilities. Height ranges from 25 to 43 inches, while electric Trendelenburg tilt offers +15° to -11° of motion. The table includes 10 inches of smooth electric longitudinal travel for precise alignment during imaging procedures.

The CFPM301 comes equipped with a multi-function hand control for simplified adjustments, a 2-inch radiolucent Comfort Foam pad, 4-inch locking casters, safety straps, and a Crescent face rest with headrest for patient comfort and stability throughout procedures.

## FEATURES

- Carbon Fiber Top - 56 in. Metal-free imaging space
- Motion 1 Electric Height Range 25" - 43" (without pad)
- Motion 2 Electric Trendelenburg tilt +15°/-11°
- Motion 3 Electric Longitudinal travel 10 in.
- Multi-function hand control
- Weight Capacity - 500 lb
- 4 in. locking casters
- Includes Patient safety straps (2)

## SPECIFICATIONS

- Patient Weight Capacity: 500 lbs.
- Carbon Fiber Top: 22" W x 84" L
- Motion 1: height range 25" - 43" (without pad)
- Motion 2: Trendelenburg tilt +15/-11°
- Motion 3: 10" longitudinal travel
- Control: 1 hand (3 motion)
- Table Weight: 475 lbs.
- Warranty: 3 years parts, 2 years labor

*\*Brochure and images are a general overview of the product. Actual product options may vary.*