

Oakworks® Fluoro Extender



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks Fluoro Extender is a radiolucent positioning accessory designed to enhance imaging accuracy during fluoroscopy procedures. Its shape accommodates various table placement applications, providing extra width where needed for optimal patient positioning. Constructed from cost-effective polycarbonate, the extender is fully radiolucent and does not interfere with imaging.

The device is simple to use and requires no attachment hardware. It slides easily under the tabletop pad, and the patient's weight holds it securely in place during the procedure. This eliminates the need for complex installation or additional fastening systems.

The Fluoro Extender is compatible with any fluoroscopy table and can be used for foot or arm support applications. Its versatile design makes it a practical optional accessory for facilities seeking to improve patient positioning and imaging outcomes.

FEATURES

- Made of a cost-effective radiolucent material - polycarbonate
- The shape is conducive to a variety of table placement applications
- Can be used with any fluoroscopy table
- Inserts under the tabletop pad and requires no attachment hardware
- For Foot or Arm Support (shown in use below)

SPECIFICATIONS

- Weight: 6 lbs
- Measurements: 23.5"W x 30" L X 1/4" Thick
- Frame: Clear Polycarbonate
- P/N: 33975

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)