

Oakworks CFPM401 Imaging Table w/Integrated Headrest



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPM401 is a 4-motion carbon fiber imaging table designed for diagnostic and interventional procedures. It features powered height adjustment (26-43 inches), Trendelenburg tilt, longitudinal travel, and lateral travel, all controlled via multi-function hand control. The table supports patients up to 500 lbs and includes an integrated headrest with crescent face rest pad.

The carbon fiber top measures 22 inches wide by 84 inches long and provides 56 inches of metal-free imaging area for unobstructed radiographic visualization. The table is equipped with a 1-inch thick radiolucent pad, two patient safety straps, and four 4-inch locking casters for mobility and stability.

The CFPM401 meets rigorous safety and regulatory standards, including FDA listing, CE marking, cETLus approval, and CB certification. It comes with a 3-year parts warranty and 2-year labor warranty.

FEATURES

- Electric 4 Movement Table: Height, Trendelenburg Tilt, Longitudinal Travel, and Lateral Travel
- Integrated Headrest
- Patient Weight Capacity: 500 lbs
- Electric Height Range: 26" to 43" (66-109 cm) without pad
- Carbon Fiber Top: 22" (56 cm) wide x 84" (214 cm) long
- 56" Metal-Free imaging space
- Hand Control (all motions)
- 1" Thick Radiolucent Table Top Pad

SPECIFICATIONS

- Electric Height Range: 26" to 44"
- Electric Trendelenburg: +15°/-11°
- Electric Longitudinal Travel: 10" (25 cm)
- Electric Lateral Travel: 8 inches (20 cm), ±4" (10 cm)
- Weight Capacity: 500 lbs (227 kg)
- Table Weight: 475 lbs (216 kg)
- Warranty: 3 years parts, 2 years labor
- Certifications: FDA listed, cETL Approved, CE Marked

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)