

Oakworks CFPMB301 Imaging Table | Bariatric Rectangular Top

FACTORY NEW

The Oakworks CFPMB301 is a bariatric-rated imaging table designed for high-quality diagnostic and interventional procedures. It features a 24" x 80" rectangular top with 56" of metal-free imaging area, supporting wheelchair transfers and full lateral C-Arm imaging mobility. The table accommodates patients up to 750 lbs.



Stock Photo - Options May Vary

The CFPMB301 offers three powered motions controlled by a single hand control: height adjustment from 26" to 44" (the industry's lowest adjustable range), Trendelenburg tilt of +15/-11°, and 10" of longitudinal travel. These capabilities enable precise positioning for a variety of imaging and procedural applications.

The table includes a 1" radiolucent pad, patient safety strap, heavy-duty 5" casters, and both foot and hand controls for all movements. It is FDA listed, cELT approved, and CE marked, with a 3-year parts and 2-year labor warranty.

FEATURES

- Rectangular Top
- Patient Weight Capacity: 750 lbs
- 3 Powered motions - Bariatric Rated
- Motion 1: Height Range 26 inches - 44 inches (66-109 cm) without pad
- Motion 2: Trendelenburg Tilt +15/-11°
- Motion 3: Longitudinal Travel 10" (25 cm)
- 56" of Metal-Free Imaging Area
- 1 Hand Control: All Motions

SPECIFICATIONS

- Patient Weight Capacity: 750 lbs
- Rectangular Bariatric Top 24" W x 80" L with 52.5" metal-free imaging space
- Motion 1: Height Range 26" - 44" (without pad)
- Motion 2: Trendelenburg Tilt: +15°/11°
- Motion 3: Lateral Tilt: ±15°
- Patient Safety Strap: 84" L with 24" of Velcro® at each end
- Table Top Pad: 1" Radiolucent Table Top Pad
- Table Weight: 475 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*