

Oakworks CFPMFXH Imaging Table with Integrated Headrest



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPMFXH is a fixed-height C-Arm imaging table designed for fluoroscopy-guided procedures in pain management and interventional settings. It features a carbon fiber top with integrated headrest, accommodating both mobile and ceiling-suspended C-Arm systems. The table provides a 56-inch metal-free imaging area for optimal image quality and C-Arm efficiency.

The table is available with either a 22" x 84" integrated headrest configuration or a 24" x 80" rectangular top, with no price difference between options. Standard height is 37", with custom heights available from 31" to 36" in one-inch increments. The 500-pound patient weight capacity makes it suitable for a wide range of clinical applications.

Mounted on a low-profile steel base with four locking casters, the table moves easily between rooms. It comes standard with hand control and includes a comfort foam pad, padded safety strap, and cranial positioning pads for patient comfort and proper positioning during procedures.

FEATURES

- Integrated Headrest
- Carbon fiber top
- Weight Capacity: 500 lbs
- Compatible with all C-arm products
- Accommodates mobile and ceiling-suspended C-Arm systems
- 56" metal-free imaging area for maximum C-Arm efficiency
- 4 locking casters
- Comes standard with hand control

SPECIFICATIONS

- Table Top Dimensions: 22" x 84" (56cm x 213cm)
- Patient weight capacity: 500 lbs (227 kg)
- Metal-free imaging area: 56" (142 cm)
- Fixed height: 37" standard
- Frame: Single Tower Steel Base
- Table weight: 415 lbs (188 kg)
- Warranty: 3 years parts, 2 years labor
- Regulatory: FDA, ETL, CE

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)