

Oakworks CFPMFXH Imaging Table - Rectangular Top



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPMFXH is a fixed-height C-Arm imaging table designed for pain management and fluoroscopy-guided procedures. It features a carbon fiber top available in either a 22" x 84" integrated headrest configuration or a 24" x 80" rectangular design, both providing a 56-inch metal-free imaging area for optimal C-Arm efficiency. The table accommodates both mobile and ceiling-suspended C-Arm systems and is compatible with all major C-Arm products.

Offered in a standard 37" height with custom height options available from 31" to 36" in one-inch increments, the table supports a 500-pound weight capacity. It is mounted on a low-profile base with four locking casters for easy mobility between rooms. The table comes standard with hand control and includes a radiolucent comfort foam pad, crescent face rest pad, supine cranial pads, and padded safety straps.

This cost-effective solution is ideal for new institutions, expanding facilities, and startup practices requiring a reliable imaging table. The fixed-height design provides stability for consistent positioning during procedures while maintaining the image quality necessary for successful interventional pain management work.

FEATURES

- Carbon fiber top
- Weight Capacity: 500 lbs (227 kg)
- Superior image quality for fluoroscopy-guided procedures
- Compatible with all major C-arm products
- Accommodates mobile and ceiling-suspended C-Arm systems
- 56" (142 cm) metal-free imaging area for maximum C-Arm efficiency
- The table is mounted on a new low-profile base
- Comes standard with hand control

SPECIFICATIONS

- Table Top Dimensions: 22" x 84" (56cm x 213cm)
- Carbon Fiber with Integrated headrest 24" x 80"
- Patient weight capacity: 500 lbs
- Metal-free imaging area: 56" (142 cm)
- Fixed height: 37" standard
- Frame: Single Tower Steel Base
- Table weight: 415 lbs (188 kg)
- Warranty: 3 years parts, 2 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*