

Oakworks Clinician Adjustable Flat Top



Stock Photo - Options May Vary

FACTORY NEW

Oakworks Clinician Adjustable Flat Top. Oakworks' Interlock leg extension system provides stability, ensuring peace of mind for both you and your client. Additionally, the manually adjustable spa table eliminates the need for electricity.

A convenient storage shelf underneath lets you store all of your massage equipment, accessories, or anything else you need close at hand. Our padding and upholstery choices ensure that your clients enjoy every second of their massage experience. Additional Features: The leg extension system uses one adjustment knob instead of two, making adjustments easier.

This adjustment process doesn't compromise the strength or integrity of the table.

FEATURES

- Padding: Plush 3"
- TerraTouch Upholstery
- Storage Shelf
- Fully Manual Table
- Width Option: 29" or 31"

SPECIFICATIONS

- Top: Flat Top
- Lift System: Adjustable
- Width: 29" or 31"
- Length: 73"
- Height Range: 24" - 36"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)