

Oakworks Clinician Adjustable Lift-Assist Backrest

FACTORY NEW

Oakworks Clinician Adjustable Lift-Assist Backrest. The Oakworks' Interlock leg extension system provides stability, ensuring peace of mind for both you and your clients.

Additionally, this manually adjustable spa table saves electricity for your facility. Moreover, the Adjustable Lift system is the most economical choice for your practice, and it includes an open shelf underneath to store all your massage equipment and accessories.

Additional Features: The Oakworks leg extension system uses one adjustment knob, making adjustments easier without compromising the table's strength or integrity. Additionally, the knob is reinforced with a nut and surrounded by a padded rubber grip that resists dirt and oil.



Stock Photo - Options May Vary

FEATURES

- Padding: Plush 3"
- TerraTouch™ Upholstery
- Storage Shelf
- Fully Manual Table
- Easy to Reach and Adjust Finger Controls even with a Patient on the Table (shown below)

SPECIFICATIONS

- Top: Backrest Top
- Lift System: Adjustable
- Width: 29" or 31"
- Length: 73"
- Height Range: 24" - 36"(with Optional 4" AeroCel™ padding)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)