

Oakworks Multi-Specialty Ultrasound Table

FACTORY NEW

The Oakworks Multi-Specialty Ultrasound Table is a versatile, ergonomically designed examination table suitable for ultrasound and sonography applications. The table features a twin tower design for strength and rigidity, with a sleek base that provides unobstructed compatibility with ultrasound equipment and promotes optimal ergonomic conditions for the sonographer.



The table offers dynamic electric height adjustment from 22 to 38 inches, supporting safe wheelchair transfers and proper working heights. It includes a bariatric-rated 550 lb lift and load capacity, along with a super comfortable 3-inch padding system available in multiple upholstery colors.

Motions include Trendelenburg and Reverse Trendelenburg positioning at ± 15 degrees, plus electric Fowler adjustment. The table is FDA Listed, cETL Approved, and CE Marked, backed by a 3-year parts warranty.

Stock Photo - Options May Vary

FEATURES

- Dynamic Electric Height Range: 22" to 38" adjustable height range for safe wheelchair transfer and ergonomic working heights
- Bariatric Rated lift and load capacity: 550 lbs.
- Twin tower design for added strength and rigidity
- The sleek base design promotes unobstructed compatibility with ultrasound equipment ensuring better ergonomic conditions for the sonographer
- Super comfortable 3" padding system
- Trendelenburg and Reverse Trendelenburg $\pm 15^\circ$
- Electric Fowler
- FDA Listed, cETL Approved, CE Marked

SPECIFICATIONS

- Patient weight capacity: 550 lbs
- Height range: 22"- 38"
- Width: 27"(option) or 30"
- Length: 72"
- Trendelenburg /Reverse: $\pm 15^\circ$
- Padding: 3" (8 cm) multi-layer
- Table weight: 359 lbs. (163 kg.)
- Certifications: FDA listed, cETL approved, CE marked

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com