

# Oakworks General Ultrasound Table w/Fowler

## FACTORY NEW

The Oakworks General Ultrasound Table with Electric Fowler is a 2-section examination table designed specifically for sonography procedures. It features an electrically adjustable Fowler section, dynamic height adjustment from 22 to 38 inches, and a 550 lb bariatric weight capacity. The sleek base design provides unobstructed compatibility with ultrasound equipment for optimal ergonomic positioning of the sonographer.



The table is constructed with a twin tower design for added strength and rigidity, and includes Trendelenburg and Reverse Trendelenburg positioning capabilities. It comes equipped with a multi-function hand control, heavy-duty locking casters, and a paper roll holder. The table features 3 inches of comfort padding and is available in 21 color options with PVC-free fabric.

Designed by leading ergonomic specialists, this table exceeds industry standards for sonography applications. It includes a 3-year parts warranty and 1-year labor warranty. The table is FDA listed, cETLus approved, CE marked, and CB certified.

*Stock Photo - Options May Vary*

## FEATURES

- Economic 2 section top table with electrically adjustable Fowler and 12 optional features to upgrade overall functionality (see Options)
- Designed by leading Ergonomic Specialists and exceeds industry standards for Sonography
- The sleek base design promotes unobstructed compatibility with ultrasound equipment ensuring better ergonomic conditions for the sonographer
- Industry-leading 3-year warranty
- Dynamic 22" to 38" adjustable height range for safe wheelchair transfer and ergonomic working heights
- 550 lb bariatric-rated lift and load capacity
- Twin tower design for added strength and rigidity
- Trendelenburg and Reverse Trendelenburg

## SPECIFICATIONS

- Patient Weight Capacity: 550 lbs. (272 kg.)
- Height Range: 22" - 38" (56-97 cm)
- Width: 27" (69 cm) or 30" (76 cm)
- Length: Fowler - 34" (86 cm), Seat - 38" (97 cm)
- Trendelenburg/Reverse Tilt:  $\pm 15^\circ$
- Padding: 3" (8 cm) multi-layer
- Table Weight: 359 lbs. (163 kg.)
- Warranty 3 Years Parts, 1 Year Labor

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)