

Oakworks Ultrasound Table Options: Side Rails

FACTORY NEW

Oakworks Ultrasound Table side rails are a set of two rails designed to enhance patient safety and stability during procedures. The rails easily tuck out of the way when not in use, providing convenient storage without taking up additional space. These side rails are a standard feature on Oakworks Vascular tables but are optional accessories for other Oakworks Ultrasound table models and must be ordered together with the table.



Stock Photo - Options May Vary

FEATURES

- A Set of 2 Side Rails
- Increase stability
- Increase security for your clients
- Easily tuck out of the way when not in use

SPECIFICATIONS

- 25" long x 6" tall (above the edge of the table)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)