

Oakworks General Ultrasound Table: 3 Section

FACTORY NEW

The Oakworks USG3 is a three-section ultrasound table designed specifically for sonography applications. It features a dynamic height range of 22" to 38" to accommodate safe wheelchair transfer and ergonomic working positions for sonographers. The table is rated for a 550 lb patient capacity and includes electric height adjustment, Trendelenburg/Reverse tilt ($\pm 15^\circ$), and electric Fowler backrest control.



Stock Photo - Options May Vary

The table's twin tower design provides structural strength and rigidity, while the sleek base design allows unobstructed compatibility with ultrasound equipment. The three-section top (Fowler 34", Seat 24", Footrest 13") can be customized with up to 12 optional features to meet specific clinical needs. A manual four-position leg rest and paper roll holder are included as standard.

The table is FDA listed, ETL approved, and CE marked. It comes with a 3-year parts warranty and 1-year labor warranty. The unit weighs 359 lbs and is available in 26" or 30" widths with a standard length of 72".

FEATURES

- Three-section top table
- 550 lb bariatric-rated lift and load capacity
- Electric Height Range 19" - 33"
- Twin tower design for added strength and rigidity
- Trendelenburg/Reverse Tilt: $\pm 15^\circ$
- Electric Fowler
- Hand Control (Height, Trendelenburg/Reverse, and Fowler)
- 3 Year Warranty

SPECIFICATIONS

- Patient weight capacity: 550 lbs.
- Height range: 22"- 38"
- Width: 26" or 30"
- Length: 72"
- Top 3 section: Fowler - 34", Seat - 24", Footrest - 13"
- Electric Trendelenburg/Reverse: $\pm 15^\circ$
- Electric Fowler Backrest: 0-75°
- Table weight: 359 lbs.

Brochure and images are a general overview of the product. Actual product options may vary.