

Oakworks Marina Table with Aqueduct Package

FACTORY NEW

Mud, Seaweed, Salt... Wet treatment tables have to be able to stand up to more than just water. The Marina is the ultimate solution for your wet table environment, accommodating a variety of heating elements and featuring the Aqueduct recessed water collection system.



Stock Photo - Options May Vary

Advanced Ergonomic Positioning for Therapist, Wet-Room Safe Battery Powered Electric Lift, and now with Full Size Fitted Wet Drape. Recessed Aqueduct: The Aqueduct is Recessed to minimize discomfort for clients and practitioners but is designed to capture almost all the water and materials used in treatments and direct them to a drain to reduce slip hazards and continuous exposure to technicians.

FEATURES

- Aqueduct™ System
- Waterproof foot control
- Charger & extra battery
- FSC® certified mahogany wood
- 3" plush padding

SPECIFICATIONS

- Top: Flat Top
- Lift System: 1 Electric Motor/Battery
- Width: 31"
- Length: 72"
- Height Range: 25" - 34"

**Brochure and images are a general overview of the product. Actual product options may vary.*