

# Oakworks Orthopedic Hi-Lo Casting Table



Stock Photo - Options May Vary

## FACTORY NEW

Oakworks Orthopedic Hi-Lo Casting Table. The Hi-Lo Casting Table is the best solution to minimize back-breaking work associated with lower extremity casting!

The 17" to 35" electric height range makes wheelchair transfer safe and easy. This table is a cost-effective component of the Ergonomic Safety Lift Plan of any medical facility. The table is FDA Listed, ETL Approved, and CE Marked.

Additional Hi-Lo Casting Table Features: 550 lb Bariatric lift capacity and load patient weight rating A unique polymer casting surface allows for easy clean-up of casting materials To view Upholstery Color Choices, see the Color Chart image in the Gallery. To customize to meet your needs, see "Options".

## FEATURES

- The best solution to minimize back-breaking work associated with lower extremity casting
- Lift-Assist 3 Section top
- 17" to 35" electric height range makes wheelchair transfer safe and easy
- A cost-effective component to any medical facility Ergonomic Safety Lift Plan
- 550 lb Bariatric lift and load patient weight rating means no more heavy lifting

## SPECIFICATIONS

- Top: 31" (79 cm) wide x 73" (185 cm) long.
- Motion 1: Electric Height 17" - 35" (43 - 89 cm).
- Motion 2: Lift-Assist Backrest 0° - 90°.
- Motion 3: Manual Leg Section 10° - 90°.
- 1 Hand Control: Controls all powered movements.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)