

# Oakworks PF Series Power Exam/Treatment Table | PF250

## FACTORY NEW

Oakworks PF Series Power Exam/Treatment Table model number: PF250. Oakworks PF Series Exam Tables give you the flexibility to create the perfect table for your needs. A wide array of options and accessories combine with great ergonomics and patient comfort to deliver great value.

The PF Series has an ADA-Compliant Height Range. Low Height Range for easy patient transfers. Oakworks PF250 Features: Powered Height, Trendelenburg (tilt), and Manual Back Rest and is designed to maximize comfort for extended treatments.

The PF250 has a Bariatric Twin Tower design has a 550 lbs. weight capacity and comes with an industry-leading 3 3-year warranty.



*Stock Photo - Options May Vary*

## FEATURES

- Motion 1: Powered Height Range
- Motion 2: Powered Tilting
- Motion 3: Manual Backrest
- Multi-function Foot Control
- Table Voltage: US (International optional)

## SPECIFICATIONS

- Motion 1: Powered 22" - 38" (56 - 97 cm) Height Range
- Motion 2: Powered +15°/-15° Tilting
- Motion 3: Manual 0-80° Backrest
- Controls: Multi-function foot control
- Table Voltage Options: US, (International optional)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*