

Oakworks PF Series Power Exam/Treatment Table | PF400

FACTORY NEW

Oakworks PF Series Power Exam/Treatment Table model number: PF400. Oakworks PF Series Exam Tables give you the flexibility to create the perfect table for your needs. A wide array of options and accessories combine with great ergonomics and patient comfort to deliver great value. The PF Series has an ADA-Compliant Height Range.



Stock Photo - Options May Vary

Low Height Range for easy patient transfers. Oakworks PF400 Features: Powered Height, Backrest, Trendelenburg (tilt), and Leg Rest, and is designed to maximize comfort for extended treatments. The PF400 has a Bariatric Twin Tower design has a 550 lbs.

weight capacity and comes with an industry-leading 3-year warranty.

FEATURES

- Motion 1: Powered Height Range
- Motion 2: Powered Tilting
- Motion 3: Powered Backrest
- Motion 4: Powered leg section
- Multi-function Foot Control

SPECIFICATIONS

- Motion 1: Powered 22" - 38" Height Range
- Motion 2: Powered +15°/-15° Tilting
- Motion 3: Powered 0-53° Backrest
- Motion 4: Powered leg section
- Controls: Multi-function Foot Control

**Brochure and images are a general overview of the product. Actual product options may vary.*