

Oakworks PG Series Power Exam/Treatment Table | PG350

FACTORY NEW

Oakworks PG Series Power Exam/Treatment Table: PG350. Oakworks, Inc. has provided ergonomic, durable equipment since 1977. It supports health and well-being globally.



Additionally, The PG Series is designed for diverse treatments with an ADA-Compliant Height Range. Oakworks PG350 Features a 3-section top that functions as a table, chair, and recliner. It includes powered height, backrest, Trendelenburg, and a manual leg rest. The table's twin tower design offers stability.

Additionally, a multi-function foot control enhances ease of use. Comfort Foam™ Padding is available in 3" or 4" options. Premium TerraTouch™ Fabric, designed by Oakworks®, is durable and long-lasting. Additionally, this PVC-free, 100% polyurethane-faced fabric features cloth backing and through-color polyurethane for enhanced durability.

Stock Photo - Options May Vary

FEATURES

- Powered Motions: Height Range, Tilting, and Backrest
- Manual Motion: 4 Position Leg Rest
- Multi-function Foot Control
- Table Voltage: US
- Patient Capacity: 550 lbs.

SPECIFICATIONS

- Motion 1: Powered 22" - 38" Height Range
- Motion 2: Powered +15°/-15° Tilting
- Motion 3: Powered 0°-75° Backrest
- Motion 4: Manual 4 position leg rest 0°-90°
- Controls: Multi-function foot control

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)