

Oakworks Physical Therapy Hi Lo Table | PT250

FACTORY NEW

Oakworks Physical Therapy Hi Lo Table: PT250. Oakworks has been providing customers with innovative medical equipment, excellent service, and cost-effective solutions for over 40 years.

Oakworks Medical has collaborated with leaders in the medical specialties we serve to create cutting-edge table designs and patient positioning devices that balance the needs of both patients and practitioners. Oakworks, Inc.

is ISO 13485 certified. The ISO standard is intended to ensure that medical device manufacturers have the required quality management systems in place to safely design, manufacture, and distribute medical devices.



Stock Photo - Options May Vary

FEATURES

- 2 Section Top
- Backrest
- Seat Section: 36 in. long
- Motion 1: Electric Height Range
- Motion 2: Manual 4 Position Head Section

SPECIFICATIONS

- Weight Capacity: 550 lb
- Product Weight: 151 lb
- 2 Section Top Dimensions: 29 in. or 31 in. wide x 72 in. long
- Flat Section: 61" long
- Head Section: 11" long

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)