

Oakworks Physical Therapy Hi Lo Table | PT300

FACTORY NEW

Oakworks Physical Therapy Hi Lo Table: PT300. Oakworks has been providing customers with innovative medical equipment, excellent service, and cost-effective solutions for over 40 years.

Oakworks Medical has collaborated with leaders in the medical specialties we serve to create cutting-edge table designs and patient positioning devices that balance the needs of both patients and practitioners. Oakworks, Inc.

is ISO 13485 certified.



Stock Photo - Options May Vary

FEATURES

- 3 Section Top: 29 in. or 31 in. wide x 73 in. long
- Backrest: 39 in. long
- Seat: 23 in. long
- Head: 11 in. long
- Motion 1: Electric Height Range

SPECIFICATIONS

- 3 Section Top: 29 in. or 31 in. wide x 73 in. long
- Backrest: 39 in. long
- Seat: 23 in. long
- Head: 11 in. long
- Motion 1: Electric Height Range 16 in. - 34 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)