

# Oakworks PowerLine™ 27" Wide Treatment Flat Top Table w/H Brace | PW27



*Stock Photo - Options May Vary*

## FACTORY NEW

The Oakworks PowerLine™ PW27 is a 27-inch wide treatment table constructed from natural bamboo, which offers superior strength and durability compared to traditional hardwood. Bamboo's tensile strength exceeds that of steel, reducing the need for frequent adjustments over the table's lifetime. The environmentally friendly material resists swelling and flaking in wet or humid conditions.

This bariatric-rated table supports patient weights up to 500 lbs and features fold-out legs for easy setup and portability. The H-brace design provides improved stability, while gusseted reinforced corners enhance durability and professional appearance. Heavy-duty vinyl upholstery resists punctures and tears, and accessory outlets at both ends offer versatility for equipment attachment.

The table arrives ready for use within minutes of setup, with assembly tools included. It measures 27 inches wide by 72 inches long with a 30-inch height and includes 2.5 inches of comfort foam padding. The PowerLine™ PW27 is backed by a full 3-year warranty and is FDA listed and CE marked.

## FEATURES

- Natural Bamboo Construction
- Bariatric Rated 500 lb. Patient Weight Capacity
- H-Brace Adding Improved Stability
- Accessory Outlets Located at both ends of the table provide versatility
- Gusseted Reinforced Corners add durability and a professional look to your table
- Heavy Duty Vinyl Resists Punctures and Tears
- Assembly Tools are Included
- Full 3-Year Warranty

## SPECIFICATIONS

- Patient Capacity: 500 lbs.
- Top Dimensions: 27"W x 72" L x 30"H
- Padding: 2.5" Comfort Foam™
- Upholstery: Premium PVC Free TerraTouch™
- H-Brace: Solid Hardwood Crossbeam Construction
- Table Weight: 125 lbs.
- Warranty: 3 years parts, 1-year labor
- P/N/ 67528

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)