

# Oakworks PowerLine™ 30" W Treatment Table w/Shelf | PW30-SHBR

**FACTORY NEW**

Oakworks PowerLine™ 30" W Treatment Table w/Shelf: PW30-SHBR.

Oakworks utilizes Natural Bamboo for its superior strength and durability under the most extreme conditions. Furthermore, environmentally friendly bamboo construction means no swelling or flaking in wet or humid conditions. Natural Bamboo has a tensile strength greater than steel and minimizes the number of adjustments typically needed over time on hardwood base tables.



*Stock Photo - Options May Vary*

**FEATURES**

- Natural Bamboo Construction
- Bariatric Rated 500 lb Patient Weight Capacity
- Full-Length Shelf Notched to Frame Adding Improved Stability
- CE Marked
- Sturdy Manual Wave Backrest Design for Incline Positioning

**SPECIFICATIONS**

- Load Weight: 500 lbs
- Height: 31"
- Width: 30"
- Length: 72"
- Top: Backrest Top

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)