

Oakworks PowerLine™ 30" Wide Treatment Table w/Backrest | PW30-BR

FACTORY NEW

The Oakworks PowerLine™ 30" Wide Treatment Table with Backrest features natural bamboo construction that provides superior strength and durability. Bamboo offers tensile strength greater than steel and resists swelling or flaking in wet or humid conditions, minimizing adjustments needed over time compared to hardwood alternatives.



Stock Photo - Options May Vary

This table is bariatric rated for 500 lb patient weight capacity and includes a sturdy manual wave backrest design for incline positioning. The fold-out legs enable easy setup, while accessory outlets at both ends provide versatility for equipment attachment. Reinforced corners enhance durability and professional appearance.

The table features thick, durable TerraTouch upholstery with superior scrape and scuff resistance, available in 21 colors. Heavy-duty vinyl resists punctures and tears, and assembly tools are included for quick setup.

FEATURES

- Natural Bamboo Construction
- Bariatric Rated 500 lb Patient Weight Capacity
- Sturdy Manual Wave Backrest Design for Incline Positioning
- Accessory Outlets Located at both ends of the table provide versatility
- Gusseted Reinforced Corners add durability and a professional look to your table

SPECIFICATIONS

- Load Weight: 500 lbs.
- Height: 30"
- Width: 30"
- Length: 72"
- Padding: 2.5" Comfort Foam™ Padding

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)