

Oakworks ProLuxe Lift-Assist: Backrest Top

FACTORY NEW

Oakworks ProLuxe Lift-Assist with Backrest Top. The ProLuxe Series is the perfect choice when you are looking for great value and excellent functionality. ProLuxe Lift-Assist Backrest Top is FSC Certified. This is an exceptional value for practitioners looking to provide simple aesthetics and bodywork, such as facials, upper body waxing, and massage.

The top is the same as found on Oakworks' most expensive tables. ProLuxe Lift-Assist Features: The manually operated backrest top provides easy adjustment for any position between flat and fully upright.



Stock Photo - Options May Vary

FEATURES

- Our largest height range
- Multi-function foot control
- Rigid twin steel beam design
- Ideal for Simple Aesthetic/bodywork such as facials, upper body waxing, and massage
- Manually Operated Backrest Top

SPECIFICATIONS

- Top: Lift-Assist Backrest Top
- Lift System: 1 Electric Motor
- Width: 29"
- Length: 73"
- Height Range: 21" - 38" (with 4" AeroCel™ padding)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)