

# Oakworks Side Lying Positioning System

## FACTORY NEW

The Oakworks Side Lying Positioning System provides ergonomic support for side-lying client positioning during massage, facials, waxing, and other therapeutic treatments. Developed with input from pregnancy instructors, the system features multiple components including a head rest, torso cushion, adjuster wedge, leg cushion, and strap that work together to ensure proper alignment and prevent neck hyper-extension.



Each component delivers targeted support for different body areas. The head rest provides cervical support and alignment, the torso cushion offers lower back and wrist support, and the leg cushion ensures proper hip support and leg relaxation. The adjuster wedge allows customization for different body types and sizes.

The system is versatile enough to function as a backrest and seated support when needed. All components feature easy-to-clean polyurethane upholstery with AeroCel foam padding, and the system includes a two-year warranty on upholstery and padding.

*Stock Photo - Options May Vary*

## FEATURES

- Helps relax and stabilize your client to make your work easier
- Stretch soft tissues and mobilize joints with more ease
- Perfect for medical massage, sports massage, and stretching
- Head Rest: Great cervical support and proper alignment support complete relaxation
- Torso Cushion: Comfort and security are ensured with this unique design
- Adjuster Wedge: The wedge creates a custom fit for maximum comfort
- Leg Cushion: Total relaxation of the leg and correct support of the hip is the goal
- Strap: Turn your flat portable table into a backrest table

## SPECIFICATIONS

- Set Includes: Head Rest, Torso Cushion, Pillow Cover (pillow not included), Adjuster Wedge, Leg Cushion, and Strap
- Upholstery: 100% polyurethane fabric, easy to clean
- Padding: AeroCel™ Foam
- P/N: 59224

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)