

Oakworks Table Option: Mounting Plate



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks mounting plate is a factory-installed accessory that enables the attachment of face rests or table extenders to compatible examination tables. For 3-section tops, the plate is available on the Fowler end only. For 2-section Fowler or rectangular ultrasound tables, it can be installed on either or both ends, with the exception of Vascular tables.

For Vascular tables equipped with Fowler, the mounting plate is available on the Fowler end. On Vascular tables without Fowler, the plate is available on the head end. Each mounting plate is ordered per end.

FEATURES

- Factory Installed
- For the 3-section top tables: Available on the Fowler end only
- For the 2-section top tables: Available on either or both ends of the tables with Fowler and rectangular ultrasound tables, except the Vascular table
- P/N: 62123

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)