

Oakworks Table Option: 12" Table Extender



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks 12-inch Table Extender lengthens compatible ultrasound tables to improve patient positioning and clinician access. It is particularly useful for carotid artery and thyroid examinations where extended reach is beneficial.

The extender is available in TerraTouch and Ultra Touch upholstery options to match existing table surfaces. It attaches to any table equipped with an Accessory Mounting Plate, which must be ordered separately.

Part number 63487.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)