

# Oakworks Ultrasound Table Options: T-Rail Adapter



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The Oakworks Ultrasound Table T-Rail Adapter enables flexible accessory mounting on ultrasound tables. The adapter can be attached at any of four designated locations on the table surface, providing multiple configuration options for different clinical setups.

This adapter is compatible with any T-Rail attachable accessory, allowing users to customize their table configuration based on specific procedural needs. The T-Rail system provides a standardized interface for quick and secure accessory installation.

The T-Rail Adapter is available as a separate purchase for existing Oakworks Ultrasound Tables, making it an economical upgrade option for facilities looking to expand their table capabilities without replacing equipment.

## **FEATURES**

- Can be attached at any of four (4) designated locations on the table
- Allows the use of any T-Rail attachable accessory
- May be ordered separately: P/N: 63336.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)