

Oakworks Urology Imaging Table: CFUR401

FACTORY NEW

The Oakworks CFUR401 is a specialized urology imaging table featuring a carbon fiber top and electronic lift tower with four powered motions: height adjustment, Trendelenburg tilt, longitudinal travel, and lateral travel. The table accommodates patients up to 500 lbs and provides a 52-inch metal-free imaging area, making it suitable for fluoroscopic procedures.



Stock Photo - Options May Vary

The table includes a removable tabletop, patient positioning straps, and 47-inch T-rails on both sides for secure patient support. It is equipped with both hand and foot controls for operator convenience, along with 4-inch locking casters for mobility and an end-mounted hand rail for safety.

Upholstery consists of medical-grade foam with TerraTouch™ fabric for durability and easy cleaning. The CFUR401 comes standard with two lateral extensions and one end extension, plus a urology bag hoop. Optional accessories such as drainage bags and stirrups are available separately.

FEATURES

- 4 Powered Motions, Height, Trendelenburg, Longitudinal travel, and Lateral travel
- Electronic Lift Tower
- Patient Weight Capacity: 500 lbs
- 52" of Metal-Free Imaging Area
- 3 Motion Controls (1 each): Foot and Hand Controls
- 47" T-rails on both sides of the top
- 4" Locking Casters
- TerraTouch Upholstery

SPECIFICATIONS

- Patient Weight Capacity: 500 lbs distributed evenly
- 26.5" to 44.5" Height Range
- Length: 47"
- +15/-11° Trendelenburg tilt
- 10" Longitudinal travel
- 52" of Metal-Free Imaging Area
- Input Service: 120VAC/15amps/60Hz
- Warranty: 3 years parts, 2 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*