

Oakworks Urology Imaging Table: CFUR301

FACTORY NEW

The Oakworks CFUR301 is an electronic urology imaging table designed for diagnostic and interventional procedures. It features three powered motions controlled by both hand and foot controls: height adjustment (26-44 inches), Trendelenburg tilt ($\pm 15/-11$ degrees), and longitudinal travel (10 inches). The table accommodates patients up to 500 lbs and includes a removable tabletop with 52 inches of metal-free imaging area.



Stock Photo - Options May Vary

The table is equipped with comprehensive positioning and safety features including patient positioning straps, 47-inch T-rails on both sides, a hand rail at the end, and 4-inch locking casters for mobility. The 1-inch medical-grade foam pad with TerraTouch fabric provides patient comfort and durability.

Additional accessories include a urology bag hoop, 2 lateral extensions, 1 end extension, and a power kill switch for safety. The electronic lift tower ensures smooth, reliable operation for clinical environments requiring precision imaging capabilities.

FEATURES

- 3 Powered Motions, Height, Trendelenburg, Longitudinal travel
- Electronic Lift Tower
- Patient Weight Capacity: 500 lbs
- 53" of Metal-Free Imaging Area
- 3 Motion Controls (1 each): Foot and Hand Controls
- Patient Positioning Straps
- 47" T-rails on both sides of the top
- 4" Locking Casters

SPECIFICATIONS

- Patient Weight Capacity: 500 lbs.
- 26" to 44" Height Range
- $+15/-11^\circ$ Trendelenburg tilt
- 10" Longitudinal travel
- Table Weight: 568 lbs
- Controls: 3 Motion Foot and 3 Motion Hand Control
- 52" of removable table metal-free imaging area
- Input Service: 120VAC/15amps/60Hz

**Brochure and images are a general overview of the product. Actual product options may vary.*