

Oakworks USWI Women's Imaging Ultrasound Table



Stock Photo - Options May Vary

FACTORY NEW

Oakworks USWI Women's Imaging Ultrasound Table. Oakworks has been providing customers with innovative medical equipment, excellent service, and cost-effective solutions for over 30 years. Furthermore, collaborating with leaders in the medical specialties resulted in cutting-edge table designs and patient positioning devices.

Moreover, Oakworks Medical is ISO 13485:2003 certified. USWI Women's Imaging Ultrasound Table Features A full-featured 3-section tabletop design for all Women's Imaging ultrasound modalities. Designed by leading Ergonomic Specialists and exceeds industry standards for Sonography. Included is an industry-leading 3-year warranty.

This ultrasound table is available in 21 beautiful colors. 22"-38" adjustable height ensures safe wheelchair transfers, ergonomic work, and supports a 550-pound bariatric lift capacity.

FEATURES

- Full-featured 3-section tabletop design for all Women's Imaging ultrasound modalities
- Designed by leading Ergonomic Specialists and exceeds industry standards for Sonography
- Industry-leading 3-year warranty
- 21 beautiful fabric colors available
- Dynamic 22" to 38" adjustable height range for safe wheelchair transfer and ergonomic working heights

SPECIFICATIONS

- 3 Section Top
- Bariatric Lift/Load Capacity: 550 lbs.
- Height Range: 22"- 38"
- Width: 27"
- Length: 72"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)