

# Oakworks Vascular Table without Fowler

## FACTORY NEW

The Oakworks Vascular Table is an electric-powered examination table designed for sonography and vascular procedures. It features a dynamic height range from 23" to 40" with Trendelenburg tilt capability (15° forward, 40° reverse), accommodating safe wheelchair transfer and ergonomic working positions. The dual tower design provides exceptional strength with a 550-pound patient weight capacity.



*Stock Photo - Options May Vary*

The table includes a split, drop-away footboard, two side rails, locking casters, and two patient restraint straps for safety and stability. All movements are controlled via a convenient hand control unit. The rectangular top measures 72" long and is available in 27" or 30" widths to suit different clinical spaces.

The table meets industry standards and regulatory requirements, with FDA listing, cETL approval, and CE marking. It comes with a comprehensive 3-year parts warranty and 1-year labor warranty, ensuring long-term reliability and support.

## FEATURES

- Electric Height Range: 23" - 40"
- Trendelenburg: 15°
- Reverse Trendelenburg: 40°
- Integrated footboard (split, drop away)
- Locking Casters
- Side Rails (2)
- Hand Control (All movements)
- 3 YEAR WARRANTY

## SPECIFICATIONS

- Patient Weight Capacity: 550 lbs
- Height Range: 23"- 40"
- Width: 27" or 30"
- Length & Top: Rectangular top - 72"
- Trendelenburg: 15°
- Reverse Trendelenburg: 40°
- Padding: 3" multi-layer
- Warranty: 3 years parts, 1-year labor

*\*Brochure and images are a general overview of the product. Actual product options may vary.*