

Pedigo Double Door OR Cabinet with Caster Base | P-8135-W/C



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Double Door Operating Room Cabinet with Caster Base model number: P-8135-W/C. Stainless steel cabinetry to maximize storage and space utilization. They are also available as free-standing or built-in models. Features 3/16" tempered glass doors with stainless steel lockable paddle latch, Double panel doors with FDA-approved polyethylene foam for glass shock resistance.

Additional Features: Adjustable nylon floor glides for leveling (freestanding models only). Adaptable to custom installations. It also features Steril-Gard construction, Radius glass panels for easy cleaning, and Adjustable stainless steel door hinges for vertical and horizontal door adjustment. It may also be field-adjusted.

Stainless steel shelves adjustable in 1/2" increments with position aids. Three-year warranty.

FEATURES

- Stainless steel cabinetry to maximize storage and space utilization
- Available as free-standing, built-in, or caster base models
- 3/16 tempered glass doors with stainless steel lockable paddle latch
- Double panel doors with FDA-approved polyethylene foam for glass shock resistance
- Adjustable nylon floor glides for leveling (free-standing models only)

SPECIFICATIONS

- Part Number: P-8135-WC
- Dimensions: 48"W x 22 1/2"D x 88 1/8"H
- Top: Sloped
- Shelves: 4
- Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)