

Pedigo Chrome Basin Stand w/ Double Stainless Steel Basins | P-79



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo Chrome Basin Stand w/Double Stainless Steel Basins: P-79 offers an ideal solution for high-efficiency environments that require multiple workspaces. Featuring a durable chrome frame, the stand ensures long-lasting strength and stability.

The open frame design allows for easy access from all sides, making it simple to perform tasks from any angle. Equipped with dual stainless steel basins, this model provides the flexibility of handling two tasks simultaneously.

The stand is mobile, with smooth-rolling casters for effortless repositioning, and it is easy to clean, offering a hygienic surface ideal for medical, dental, or laboratory applications.

FEATURES

- Chrome Stand
- Two Stainless Steel Basins
- Open frame design for easy access from any side
- Double ball bearing swivel casters
- Mobile

SPECIFICATIONS

- Overall Dimensions: 14 D x 28"W x 34"H
- Two removable 13¼" diameter, 7-quart stainless steel basins
- Stand Construction: Chrome
- Basin Construction: Stainless Steel
- Casters: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)