

Pedigo Chrome Basin Stand w/Single Stainless Steel Basin | P-78



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo Chrome Basin Stand w/Single Basin: P-78 features a durable chrome frame with an open design for easy access from any side. The mobile stand is equipped with smooth casters for easy movement and is designed for easy cleaning and maintenance.

The stainless steel basin offers a hygienic and durable surface, making it ideal for medical, dental, and laboratory environments. Its sleek, functional design ensures both convenience and reliability.

FEATURES

- Chrome Stand
- Stainless Steel Basin
- Open frame design for easy access from any side
- Available in single and double-basin designs
- 2" diameter double ball bearing swivel casters

SPECIFICATIONS

- Overall Dimensions: 14 D x 14"W x 34"H
- Basin diameter: One removable 13¼"
- Basin Capacity: 7 quarts
- Stand Construction: Chrome
- Basin Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)