

Pedigo Chrome/All Welded Steel Footstool | P-10



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Chrome/All-Welded Steel Footstool: P-10. The Pedigo P-10 offers strength and stability with its all-welded steel construction. The durable 16-gauge tubular steel legs provide lasting support.

A corrugated rubber top surface provides secure footing, while rubber floor tips protect your floors and prevent slipping. Pedigo footstools are tip-resistant for added safety and are easy to maintain, making them a reliable and practical choice for various settings, including healthcare environments and offices.

FEATURES

- Strong, all-welded steel
- Legs constructed of 1 O.D. 16-gauge tubular steel
- Tip resistant, rugged, and easy to maintain
- Corrugated black rubber top surface
- Rubber floor tips

SPECIFICATIONS

- Construction: All-welded Steel
- Dimensions: 12"D x 16"W x 8¼"H
- Weight Capacity: 300 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)