

Pedigo Closed Case Cart w/Complete Accessories | CDS-242-C



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Closed Case Cart w/Complete Accessories: CDS-242-C. Pedigo's Surgical Case Carts are constructed from a corrosion-resistant grade of stainless steel and built to STERIL-GARD Specifications. STERIL-GARD is the concept and process of manufacturing products to minimize or eliminate any areas where dirt, fluid, and bacterial debris can become entrapped.

Additional Features: Designed to transport sterile packs from Central Supply to Surgery. You can modify these carts to fit elevators and automatic cart pick-up systems.

The carts feature fully-welded stainless steel construction, sound-proofed tops and doors for quiet operation, and exterior push handles for better control and more interior space.

FEATURES

- All carts built to Steril-Gard specifications
- Fully welded all stainless steel construction
- Sound-proofed tops and doors for quiet operation
- Exterior push handles for better control and more interior space
- One-piece non-marring continuous wrap-around bumper

SPECIFICATIONS

- Model number CDS-242-C
- Cabinet Size: 40¼" W x 24¼" x 30"H
- Overall Dimensions: 45 7/8" W x 28 5/8" D x 39"H
- Construction: Stainless Steel
- Doors: Sealed

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)