

Pedigo Closed Case Cart w/Complete Accessories | CDS-245-C



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Closed Case Cart w/Complete Accessories: CDS-245-C. Pedigo's Surgical Case Carts are constructed from a corrosion-resistant grade of stainless steel and built to STERIL-GARD Specifications.

STERIL-GARD is the concept and process of manufacturing products to minimize or eliminate any areas where dirt, fluid, and bacterial debris can become entrapped. Additional Features: Designed to transport sterile packs from Central Supply to Surgery.

You can modify these carts to fit elevators and automatic cart pick-up systems. The carts feature fully-welded stainless steel construction, sound-proofed tops and doors for quiet operation, and exterior push handles for better control and more interior space.

FEATURES

- All carts built to Steril-Gard specifications
- Fully welded all stainless steel construction
- Sound-proofed tops and doors for quiet operation
- Exterior push handles for better control and more interior space
- One-piece non-marring continuous wrap-around bumper
- Positive lock door latches

SPECIFICATIONS

- Model number: CDS-245-C
- Cabinet Size: 40 $\frac{1}{4}$ "W x 24 $\frac{1}{4}$ " x 47"H
- Overall Dimensions: 45 $\frac{7}{8}$ "W x 28 $\frac{5}{8}$ "D x 56"H
- Construction: Stainless Steel
- Doors: Sealed

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)