

# Pedigo Deluxe Lab Stool | P-55

## FACTORY NEW

Pedigo Deluxe Lab Stool model number P-55. Today, as Pedigo celebrates 66 years of service and three generations of family ownership, they continue to lead the way in product innovation and increased sterility control.

Pedigo is committed to superior customer service. Custom upholstery colors are available to match every decor.

Large inventory for fast delivery. Pedigo P-55 Features: Designed for easy cleaning, the unit has screw shaft height adjustment, plastic glide, and an upholstered spring back.



Stock Photo - Options May Vary

## FEATURES

- Deluxe Lab Stool
- Screw shaft height adjustment
- Plastic Glide
- Upholstered spring back
- Cushioned seat

## SPECIFICATIONS

- Screw shaft height adjustment: 21-5/8" to 33-5/8"
- Square Cushioned Seat Dimensions: 14" x 15"
- Weight Capacity: 300 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)