

Pedigo Double Pod ivNow® Modular Thermal Fluid Warmer | P-2302-A



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-2302-A Double Pod ivNow® Thermal Fluid Warmer. Pedigo's ivNow modular fluid warmers provide warmed fluid at your fingertips, exactly where and when you need it. These small, convenient units offer flexibility, allowing placement in various combinations and environments.

Countertop, wall-mounted, and mobile units are all available depending on your requirements. Features Include: ivNow quickly warms and maintains safe intravenous and irrigation fluids temperatures.

Each specially contoured warming surface cradles one (1) liter solution bag. Each cavity has its own control and LED display, allowing for precise and independent temperature management.

FEATURES

- Pedigo's ivNow fluid warmers are small, convenient units
- Can be placed in a number of different combinations and environments
- Countertop, wall mounted and mobile units are all available depending on your particular requirements
- Quickly warms and maintains safe temperatures of intravenous and irrigation fluids
- Each specially contoured warming surface cradles one (1) liter solution bag

SPECIFICATIONS

- Dimensions: 13.4" H x 14.5"W
- Capacity: 2 bags (liter)
- Display: LED
- Temp Sensors: 2
- Accuracy: Within $\pm 0.2^{\circ}\text{C}$ ($\pm 0.3^{\circ}\text{F}$) of the setpoint temperature

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)