

# Pedigo Elite Series Warming Cabinet | P-2210



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo P-2210 Elite Series Warming Cabinet. Pedigo leads in designing and constructing high-quality stainless steel and chrome products. Their products are used in emergency rooms, exam rooms, operating rooms, and sterile processing departments..

And after more than a quarter of a century, the Pedigo name continues to be synonymous with quality, durability, sterility, and service to the customer. Pedigo Elite Series Warming Cabinets Feature: Single-chamber warming cabinet with White epoxy-coated steel interior insert and exterior casing.

Single pane energy-efficient e-coated glass window door allows inventory observation and easy, hands-free door release. The door has a full gasket and is hinged on the right side of the unit.

## FEATURES

- Single-chamber warming cabinet
- White epoxy-coated steel interior insert and exterior casing
- Single pane energy-efficient e-coated glass window door allows inventory observation.
- Easy, hands-free door release
- The door has a full gasket and hinged on the right side of the unit

## SPECIFICATIONS

- Temperature Range: 98°-180°
- Compartment Capacity: 2.5 cubic feet
- Compartment Dimensions: 14½" W x 18½" D x 14-5/8" H
- Overall Dimensions: 18½" W x 22-13/16"D x 21-7/8" H
- Single-chamber warming cabinet

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)