

Pedigo Elite Series Blanket Warmer | P-2220



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-2220 Elite Series Blanket Warmer. The Pedigo P-2220 Elite Series Blanket Warmer delivers dependable performance for busy medical environments. This single-chamber warmer heats blankets quickly and evenly. It uses an efficient, balanced heating system that lowers energy use and reduces heat loss.

Clinicians gain reliable warmth for patient comfort, and staff can manage temperatures with confidence. Furthermore, the cabinet supports Fahrenheit and Celsius settings for flexible operation. You control temperatures from 98° to 180°F, which ensures ideal warmth for various clinical needs.

The cabinet also includes a dedicated warming shut-off system that prevents overheating without relying on electronic controls. As a result, your team maintains safety during every shift.

FEATURES

- Electronic control
- Efficient, the balanced heating system allows for low energy consumption and minimal heat loss
- Adjustable temperature range of 98 - 180°F (37 - 82°C)
- Operate in Fahrenheit or Celsius
- Four-digit L.E.D. display

SPECIFICATIONS

- Compartment Capacity: 3.5 cubic feet
- Compartment Dimensions: 14½"W x 18½"D x 19½"H
- Overall Dimensions: 18½"W x 20-5/8"□D x 27-7/8"H
- Single-chamber warming cabinet
- Construction: White epoxy-coated steel interior insert and exterior casing

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)