

Pedigo Free-Standing Double Door OR Cabinet | P-8140



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Free-Standing Double Door Operating Room Cabinet model number: P-8140. Stainless steel cabinetry to maximize storage and space utilization. 3/16" tempered glass doors with stainless steel lockable paddle latch and Double panel doors with FDA-approved polyethylene foam for glass shock resistance. Additional Features: Adjustable nylon floor glides for leveling (freestanding models only).

Adaptable to custom installations, Steril-Gard construction, Radius glass panels for easy cleaning, and Adjustable stainless steel door hinges for vertical and horizontal door adjustment. They are adjusted in the field. Stainless steel shelves are adjustable in 1/2" increments with position aids—three-year warranty.

Note: to view additional carts, click the following link [Case Carts & Operating Room Cabinets](#).

FEATURES

- Stainless steel cabinetry to maximize storage and space utilization.
- Available with caster base or built-in
- 3/16" tempered glass doors with stainless steel lockable paddle latch
- Double panel doors with FDA-approved polyethylene foam for glass shock resistance
- Adjustable nylon floor glides for leveling (freestanding models only)

SPECIFICATIONS

- Part Number: P-8140
- Overall Dimensions: 47" x 24" x 72"
- Top: Flat
- Shelves: 4
- Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*