

Pedigo Hydraulic Surgeon's Stool | P-6000



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Hydraulic Surgeon's Stool: P-6000. Hydraulic Surgeon Chairs/Stools are ideal for surgery, anesthesiology, ENT, and reconstructive procedures.

STERIL-GARD is the concept and process of manufacturing products to minimize or eliminate any areas where dirt, fluid, and bacterial debris can become entrapped.

Pedigo P-6000 Additional Features: All Pedigo products are easier to clean, longer-lasting, and, most importantly - reduce the transmission of bacteria. The P-6000 features include: foot-operated hydraulic height adjustment, Single pedal height adjustment from 22-3/4" to 30-1/4", cushioned 16" seat, 21" diameter base, and Powder-coated five caster base. In addition, it comes complete with high-grade 3" casters - 3 casters are locking.

FEATURES

- Foot-Operated Hydraulic Height Adjustment
- Single Pedal Height Adjustment from 22-3/4" to 30-1/4"
- 21" Diameter Base
- Powder-Coated 5 Caster Base
- High Grade 3" Casters

SPECIFICATIONS

- Base Diameter: 21"
- Seat Diameter: 16"
- Seat Height Adjustment Range: 22-3/4" to 30-1/4"
- Casters Size: 3"
- Locking Casters: 3

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)